


Aquatrainer



HORSE GYM

► We create champions ◀ **2000**



www.horse-gym-2000.de

Basic principle and fields of application

The Aquatrainer is the latest product of the successful range of HORSE-GYM-2000 treadmills for horses. It allows to make horses move within the scope of a therapy without charging them, while observing their movements. Movement in the water is unexcelled with regard to the training of muscles and fitness. Indeed, the force of gravity is almost neutralized by the water and the muscles and joints thus are not loaded with the full weight. Water offers a soft resistance to increase the power. Such resistance is increased by faster movements, so that the natural load is further increased. The Aquatrainer can be filled with water individually for each horse up to a height of 1.20 m. The Aquatrainer is used for rehabilitation purposes and for supporting the training of competition/race horses and horses for leisure riding.



AQUA-TRAINER

Application and advantages

- Rehabilitation after surgeries
- Rehabilitation after sinew and ligament problems
- Treatment of spinal and back disorders
- Rehabilitation of performance capacity
- General revitalization
- Supplementary rehabilitation training
- Release of muscles and joints
- Cooling effect/stimulation

- Massage of muscles
- Healthy blood flow
- Build-up of muscles

plus:

- Simple maintenance (among others, the belt can be completely removed from the reservoir)
- Low water consumption
- Easy cleaning





HORSE-GYM 2000

Quality and safety

HORSE-GYM-2000 products excell from other products especially with regard to quality and safety, since all details are of high-quality only.

For example, a light barrier stops the belt when the horse falls and immediately empties the water reservoir.

Equally, the horse's easy mounting on and descending from the belt as well as the step-by-step starting of the belt and filling of the reservoir allows you to work without stress. The safety of human beings and animals is of utmost importance to us.



WELLNESS FOR YOU

Function and process

The tried and tested HORSE-GYM-2000 technology of the belt floor is also used in the Aquatrainer. It allows the horse to equally and planely step on the floor the quality of which is similar to a naturally grown grass or forest floor.

The adjustment of the belt speed and the water level in the reservoir allows a well-aimed adjustment to the animal's requirements. Depending on the requirements and the horse's condition, the Aquatrainer is filled with water by means of a pumping system and the horse is running in the water.



Relaxation, gymnastification and strengthening of the back and leg muscles with soft movements and cooling



Water level can be increased by 10 cm increments up to 1.20 m at the push of a button.



To be used even without water as a standard treadmill offering all benefits.

R PRECIOUS ANIMALS

Technology and filtering

Only single components meeting maximum requirements that were tested in detail are used for our treadmills and water treadmills. The same is true for the Aquatrainer's unique filtering system.

The Aquatrainer's stimulating effect is also shown by the horse's defecation. Thanks to an elaborated, low-maintenance filter system, such impurities are removed immediately. Thus the water can be reused again and again. For special hygiene purposes, such as for clinical applications, we offer special disinfection systems.



Impurities such as horse dung are removed thanks to a self-cleaning, low-maintenance permanent filter.



The water is reused.



The water's quality and level are always visible.

HORSE-GYM-2000 TREA

Your individual Aquatrainer

Each Aquatrainer made by our company is customized to your individual needs. Thus, we are able to respond to your local needs.

The filter system can be connected according to the space available, at the left or right side of the reservoir, without extra charge. The filter system could also be installed in a different room than the Aquatrainer reservoir.

Various basic models are available:

Aquatrainer model W2:

- small reservoir with a water level of up to 1 m
- 2 doors
(for even simpler mounting and descending of the horses)

Aquatrainer model W3:


- big reservoir with a water level of up to 1.20 m
(the horse gets even more lift: back and legs are thus even more released)
- 2 doors

Standard:

- L x W x H: 600 cm x 230 cm x 200 cm
- Weight: 2000 kg
- Power supply: 5 kW / 380 V
- Water tank 1 x 4000 l
- Filtering system for suspended matter
- to be connected on both sides
- Standard piping
- Control flaps
- Small components

Options:

- individual painting
- cooling / heating unit
- 6000 L tank
- additional disinfection upon request
- salt water coating
- single phase





HORSE GYM
▶ We create champions ◀ **2000**

Address:

HORSE-GYM-2000 GmbH

Owner: Siegfried Mitzel

Langgasse 18

86655 Harburg-Großsorheim

Germany

WE MAKE CHAMPIONS – RECENTLY EVEN WORLD CHAMPIONS!


Visit us on site:

Großsorheim is an incorporation of the town of Harburg, located in South Germany only 1.5 hours drive away from Munich, Stuttgart or Ulm. Thanks to its vicinity to the A6, A7, A8, A9 motorway, you'll quickly reach us no matter from what direction you come! Contact us today to make an individual appointment with us!

We look forward to consulting you!

Tel.: +49 (0) 90 80 / 92 22 99

Fax: +49 (0) 90 80 / 92 26 32



You'll find more information at

E-Mail: info@horse-gym-2000.de

www.horse-gym-2000.de